

# Thioredoxin-dependent peroxide reductase, mitochondrial

Catalog No: #AP87838

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87838-1 50ug #AP87838-2 100ug #AP87838-3 1mg

## Description

Product Name	Thioredoxin-dependent peroxide reductase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Antioxidant protein 1,AOP-1,PRX III,Perioredoxin-3,Protein MER5,Aop1,Mer5
Accession No.	Uniprot ID:P20108
Uniprot	P20108
GeneID	11757;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Prdx3

Note: This product is for in vitro research use only