

Proteoglycan 4

Catalog No: #AP87839



Package Size: #AP87839-1 50ug #AP87839-2 100ug #AP87839-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteoglycan 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lubricin, Megakaryocyte-stimulating factor, Superficial zone proteoglycan, Msf, Szp
Accession No.	Uniprot ID:Q9JM99
Uniprot	Q9JM99
GeneID	96875;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Prg4

Note: This product is for in vitro research use only