

Semaphorin-6B

Catalog No: #AP87840



Package Size: #AP87840-1 50ug #AP87840-2 100ug #AP87840-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Semaphorin-6B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Semaphorin-Z, Sema Z, SEMAN, SEMAZ, UNQ1907/PRO4353
Accession No.	Uniprot ID:Q9H3T3
Uniprot	Q9H3T3
GeneID	10501;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SEMA6B

Note: This product is for in vitro research use only