

D-3-phosphoglycerate dehydrogenase

Catalog No: #AP87841



Package Size: #AP87841-1 50ug #AP87841-2 100ug #AP87841-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	D-3-phosphoglycerate dehydrogenase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	3-PGDH,2-oxoglutarate reductase,Malate dehydrogenase
Accession No.	Uniprot ID:A5GFY8
Uniprot	A5GFY8
GeneID	100144529;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PHGDH

Note: This product is for in vitro research use only