

Selenoprotein S

Catalog No: #AP87844



Package Size: #AP87844-1 50ug #AP87844-2 100ug #AP87844-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenoprotein S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	SelS,SELS
Accession No.	Uniprot ID:F1SR90
Uniprot	F1SR90
GeneID	100151836;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SELENOS

Note: This product is for in vitro research use only