

Transmembrane 4 L6 family member 20

Catalog No: #AP87848



Package Size: #AP87848-1 50ug #AP87848-2 100ug #AP87848-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Transmembrane 4 L6 family member 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q3T0Z4
Uniprot	Q3T0Z4
GeneID	767902;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TM4SF20

Note: This product is for in vitro research use only