

Secreted frizzled-related protein 5

Catalog No: #AP87850



Package Size: #AP87850-1 50ug #AP87850-2 100ug #AP87850-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Secreted frizzled-related protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	sFRP-5,Frizzled-related protein 1b,FRP-1b,Secreted apoptosis-related protein 3,SARP-3,FRP1B,SARP3
Accession No.	Uniprot ID:Q5T4F7
Uniprot	Q5T4F7
GeneID	6425;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SFRP5

Note: This product is for in vitro research use only