

Sequestosome-1

Catalog No: #AP87852



Package Size: #AP87852-1 50ug #AP87852-2 100ug #AP87852-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sequestosome-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Protein kinase C-zeta-interacting protein,PKC-zeta-interacting protein,Ubiquitin-binding protein p62,Zip
Accession No.	Uniprot ID:O08623
Uniprot	O08623
GeneID	113894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sqstm1

Note: This product is for in vitro research use only