

Sperm-associated antigen 1

Catalog No: #AP87853



Package Size: #AP87853-1 50ug #AP87853-2 100ug #AP87853-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sperm-associated antigen 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HSD-3.8, Infertility-related sperm protein Spag-1
Accession No.	Uniprot ID:Q07617
Uniprot	Q07617
GeneID	6674;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SPAG1

Note: This product is for in vitro research use only