

Runt-related transcription factor 1

Catalog No: #AP87854

Package Size: #AP87854-1 50ug #AP87854-2 100ug #AP87854-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Runt-related transcription factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Acute myeloid leukemia 1 protein,Core-binding factor subunit alpha-2,CBF-alpha-2,Oncogene AML-1,Polyomavirus enhancer-binding protein 2 alpha B subunit,PEA2-alpha B,SL3-3 enhancer factor 1 alpha B subunit,SL3/AKV core-binding factor alpha B subunit,Aml1,Cbfa2,Pebp2ab
Accession No.	Uniprot ID:Q03347
Uniprot	Q03347
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Runx1

Note: This product is for in vitro research use only