

Retinoic acid receptor RXR-gamma

Catalog No: #AP87855



Package Size: #AP87855-1 50ug #AP87855-2 100ug #AP87855-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor RXR-gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Nuclear receptor subfamily 2 group B member 3, Retinoid X receptor gamma, NR2B3
Accession No.	Uniprot ID:Q0GFF6
Uniprot	Q0GFF6
GeneID	445465;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RXRG

Note: This product is for in vitro research use only