

Mothers against decapentaplegic homolog 7

Catalog No: #AP87856



Package Size: #AP87856-1 50ug #AP87856-2 100ug #AP87856-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAD homolog 7, Mothers against decapentaplegic homolog 8, MAD homolog 8, SMAD family member 7, SMAD 7, MADH7, MADH8
Accession No.	Uniprot ID:O15105
Uniprot	O15105
GeneID	4092;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SMAD7

Note: This product is for in vitro research use only