

Serologically defined colon cancer antigen 3

Catalog No: #AP87859



Package Size: #AP87859-1 50ug #AP87859-2 100ug #AP87859-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serologically defined colon cancer antigen 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Antigen NY-CO-3
Accession No.	Uniprot ID:Q96C92
Uniprot	Q96C92
GeneID	10807;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SDCCAG3

Note: This product is for in vitro research use only