

## CD2L5 Antibody

Catalog No: #44458



Package Size: #44458-1 50ul #44458-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD2L5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CD2L5 Antibody detects endogenous levels of total CD2L5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CD2L5
Target Name	CD2L5
Other Names	CDK13, CDC2L, CDC2L5, CHED, Cyclin-dependent kinase 13, HCDK13, KIAA1791, CDC2-related protein kinase 5
Accession No.	Swiss-Prot#: Q14004NCBI Gene ID:8621
Uniprot	Q14004
GeneID	8621;
Calculated MW	165 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only