

Selenium-binding protein 1

Catalog No: #AP87864



Package Size: #AP87864-1 50ug #AP87864-2 100ug #AP87864-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Selenium-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	56 kDa selenium-binding protein, SBP56, Selenium-binding protein 2, Sbp
Accession No.	Uniprot ID:Q8VIF7
Uniprot	Q8VIF7
GeneID	103689947;140927;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Selenbp1

Note: This product is for in vitro research use only