

Estrogen sulfotransferase

Catalog No: #AP87867



Package Size: #AP87867-1 50ug #AP87867-2 100ug #AP87867-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Estrogen sulfotransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Sulfotransferase 1E1,ST1E1,Sulfotransferase, estrogen-preferring,OST,STE
Accession No.	Uniprot ID:P19217
Uniprot	P19217
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SULT1E1

Note: This product is for in vitro research use only