

Suppressor of cytokine signaling 1

Catalog No: #AP87870



Package Size: #AP87870-1 50ug #AP87870-2 100ug #AP87870-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Suppressor of cytokine signaling 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SOCS-1,JAK-binding protein,JAB,STAT-induced STAT inhibitor 1,SSI-1,Tec-interacting protein 3,TIP-3,SSI1,TIP3
Accession No.	Uniprot ID:O15524
Uniprot	O15524
GeneID	8651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SOCS1

Note: This product is for in vitro research use only