

Transcription factor SOX-18

Catalog No: #AP87874



Package Size: #AP87874-1 50ug #AP87874-2 100ug #AP87874-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor SOX-18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Sox-18
Accession No.	Uniprot ID:P43680
Uniprot	P43680
GeneID	20672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Sox18

Note: This product is for in vitro research use only