

Signal recognition particle 54 kDa protein

Catalog No: #AP87875



Package Size: #AP87875-1 50ug #AP87875-2 100ug #AP87875-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle 54 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	SRP54
Accession No.	Uniprot ID:P61010
Uniprot	P61010
GeneID	403953;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRP54

Note: This product is for in vitro research use only