

Sperm-specific antigen 2

Catalog No: #AP87876



Package Size: #AP87876-1 50ug #AP87876-2 100ug #AP87876-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sperm-specific antigen 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cleavage signal-1 protein,CS-1,Ki-ras-induced actin-interacting protein,CS1,KIAA1927,KRAP
Accession No.	Uniprot ID:P28290
Uniprot	P28290
GeneID	6744;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SSFA2

Note: This product is for in vitro research use only