

Somatostatin receptor type 2

Catalog No: #AP87877

Package Size: #AP87877-1 50ug #AP87877-2 100ug #AP87877-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SS-2-R,SRIF-1
Accession No.	Uniprot ID:P30874
Uniprot	P30874
GeneID	6752;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SSTR2

Note: This product is for in vitro research use only