

## Estrogen sulfotransferase

Catalog No: #AP87878



Package Size: #AP87878-1 50ug #AP87878-2 100ug #AP87878-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Estrogen sulfotransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	EST-1,Sulfotransferase 1E1,ST1E1,Sulfotransferase, estrogen-preferring,STE
Accession No.	Uniprot ID:P49888
Uniprot	P49888
GeneID	6783;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SULT1E1

Note: This product is for in vitro research use only