

Signal transducer and activator of transcription 4

Catalog No: #AP87879



Package Size: #AP87879-1 50ug #AP87879-2 100ug #AP87879-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q14765
Uniprot	Q14765
GeneID	6775;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STAT4

Note: This product is for in vitro research use only