

Serine/arginine-rich splicing factor 2

Catalog No: #AP87882



Package Size: #AP87882-1 50ug #AP87882-2 100ug #AP87882-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/arginine-rich splicing factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Splicing component, 35 kDa, Splicing factor SC35, SC-35, Splicing factor, arginine/serine-rich 2, Srsf2
Accession No.	Uniprot ID:Q6PDU1
Uniprot	Q6PDU1
GeneID	494445;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Srsf2

Note: This product is for in vitro research use only