

Sulfotransferase 1A1

Catalog No: #AP87886



Package Size: #AP87886-1 50ug #AP87886-2 100ug #AP87886-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sulfotransferase 1A1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	ST1A1,Aryl sulfotransferase,Aryl sulfotransferase IV,ASTIV,Minoxidil sulfotransferase,Mx-ST,PST-1,Phenol sulfotransferase,Sulfokinase,Tyrosine-ester sulfotransferase,St1a1
Accession No.	Uniprot ID:P17988
Uniprot	P17988
GeneID	83783;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sult1a1

Note: This product is for in vitro research use only