

Kit ligand

Catalog No: #AP87887



Package Size: #AP87887-1 50ug #AP87887-2 100ug #AP87887-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Kit ligand |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Chicken |
| Other Names | Mast cell growth factor, MGF, Stem cell factor, SCF, c-Kit ligand, SCF |
| Accession No. | Uniprot ID:Q09108 |
| Uniprot | Q09108 |
| GeneID | 396028; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: KITLG

Note: This product is for in vitro research use only