

Signal recognition particle 54 kDa protein

Catalog No: #AP87889



Package Size: #AP87889-1 50ug #AP87889-2 100ug #AP87889-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle 54 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SRP54
Accession No.	Uniprot ID:Q6AYB5
Uniprot	Q6AYB5
GeneID	116650;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Srp54

Note: This product is for in vitro research use only