

## Ribonucleoside-diphosphate reductase large subunit

Catalog No: #AP87892



Package Size: #AP87892-1 50ug #AP87892-2 100ug #AP87892-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ribonucleoside-diphosphate reductase large subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ribonucleoside-diphosphate reductase subunit M1,Ribonucleotide reductase large subunit
Accession No.	Uniprot ID:P07742
Uniprot	P07742
GeneID	20133;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Rrm1

Note: This product is for in vitro research use only