

Sphingosine 1-phosphate receptor 1

Catalog No: #AP87893



Package Size: #AP87893-1 50ug #AP87893-2 100ug #AP87893-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Sphingosine 1-phosphate receptor 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | S1P receptor 1,Endothelial differentiation G-protein coupled receptor 1,Sphingosine 1-phosphate receptor Edg-1,S1P receptor Edg-1,CD363,CHEG1,EDG1 |
| Accession No. | Uniprot ID:P21453 |
| Uniprot | P21453 |
| GeneID | 1901; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S1PR1

Note: This product is for in vitro research use only