

Sodium/glucose cotransporter 2

Catalog No: #AP87894



Package Size: #AP87894-1 50ug #AP87894-2 100ug #AP87894-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sodium/glucose cotransporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Na(+)/glucose cotransporter 2,Low affinity sodium-glucose cotransporter,Solute carrier family 5 member 2,Sgt2
Accession No.	Uniprot ID:P53792
Uniprot	P53792
GeneID	64522;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc5a2

Note: This product is for in vitro research use only