

SH3 domain-binding glutamic acid-rich-like protein 3

Catalog No: #AP87895



Package Size: #AP87895-1 50ug #AP87895-2 100ug #AP87895-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SH3 domain-binding glutamic acid-rich-like protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SH3 domain-binding protein 1,SH3BP-1,P1725
Accession No.	Uniprot ID:Q9H299
Uniprot	Q9H299
GeneID	83442;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SH3BGRL3

Note: This product is for in vitro research use only