

NAD-dependent protein deacetylase sirtuin-7

Catalog No: #AP87896



Package Size: #AP87896-1 50ug #AP87896-2 100ug #AP87896-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent protein deacetylase sirtuin-7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Regulatory protein SIR2 homolog 7,SIR2-like protein 7,Sir2l7
Accession No.	Uniprot ID:Q8BKJ9
Uniprot	Q8BKJ9
GeneID	209011;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sirt7

Note: This product is for in vitro research use only