

# Endophilin-A2

Catalog No: #AP87898



Package Size: #AP87898-1 50ug #AP87898-2 100ug #AP87898-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endophilin-A2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endophilin-2,SH3 domain protein 2B,SH3 domain-containing GRB2-like protein 1,SH3p8,Sh3d2b
Accession No.	Uniprot ID:Q62419
Uniprot	Q62419
GeneID	20405;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sh3gl1

Note: This product is for in vitro research use only