

# NAD-dependent protein deacylase sirtuin-5, mitochondrial

Catalog No: #AP87900

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87900-1 50ug #AP87900-2 100ug #AP87900-3 1mg

## Description

|                    |  |
|--------------------|--|
| Product Name       | NAD-dependent protein deacylase sirtuin-5, mitochondrial |
| Brief Description  | Recombinant Protein                                      |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                             |
| Species Reactivity | Rt   |
| Other Names        | Regulatory protein SIR2 homolog 5,SIR2-like protein 5    |
| Accession No.      | Uniprot ID:Q68FX9  |
| Uniprot            | Q68FX9   |
| GeneID             | 306840;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Sirt5

Note: This product is for in vitro research use only