

40S ribosomal protein S14

Catalog No: #AP87902



Package Size: #AP87902-1 50ug #AP87902-2 100ug #AP87902-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	40S ribosomal protein S14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Small ribosomal subunit protein uS11,PRO2640
Accession No.	Uniprot ID:P62263
Uniprot	P62263
GeneID	6208;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RPS14

Note: This product is for in vitro research use only