

Serine/threonine-protein kinase SIK1

Catalog No: #AP87903



Package Size: #AP87903-1 50ug #AP87903-2 100ug #AP87903-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase SIK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Salt-inducible kinase 1, SIK-1, Serine/threonine-protein kinase SNF1-like kinase 1, Serine/threonine-protein kinase SNF1LK, SIK, SNF1LK
Accession No.	Uniprot ID:P57059
Uniprot	P57059
GeneID	102724428;150094;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SIK1

Note: This product is for in vitro research use only