

Radical S-adenosyl methionine domain-containing protein 2

Catalog No: #AP87904

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87904-1 50ug #AP87904-2 100ug #AP87904-3 1mg

Description

Product Name	Radical S-adenosyl methionine domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytomegalovirus-induced gene 5 protein,Viperin,Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible,CIG5
Accession No.	Uniprot ID:Q8WXG1
Uniprot	Q8WXG1
GeneID	91543;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RSAD2

Note: This product is for in vitro research use only