

Antigen peptide transporter 2

Catalog No: #AP87906



Package Size: #AP87906-1 50ug #AP87906-2 100ug #AP87906-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Antigen peptide transporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APT2,ATP-binding cassette sub-family B member 3,Peptide supply factor 2,Peptide transporter PSF2,PSF-2,Peptide transporter TAP2,Peptide transporter involved in antigen processing 2,Really interesting new gene 11 protein,ABCB3,PSF2,RING11,Y1
Accession No.	Uniprot ID:Q03519
Uniprot	Q03519
GeneID	6891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TAP2

Note: This product is for in vitro research use only