

Bile salt sulfotransferase 1

Catalog No: #AP87907



Package Size: #AP87907-1 50ug #AP87907-2 100ug #AP87907-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bile salt sulfotransferase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hydroxysteroid sulfotransferase,ST,Sulfotransferase 2A1,ST2A1,Sta1
Accession No.	Uniprot ID:P52843
Uniprot	P52843
GeneID	20859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sult2a1

Note: This product is for in vitro research use only