

Extracellular sulfatase Sulf-1

Catalog No: #AP87908



Package Size: #AP87908-1 50ug #AP87908-2 100ug #AP87908-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Extracellular sulfatase Sulf-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RSulfFP1,Sulfatase FP
Accession No.	Uniprot ID:Q8VI60
Uniprot	Q8VI60
GeneID	171396;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Sulf1

Note: This product is for in vitro research use only