

Sestrin-1

Catalog No: #AP87910



Package Size: #AP87910-1 50ug #AP87910-2 100ug #AP87910-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sestrin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p53-regulated protein PA26,PA26,SEST1
Accession No.	Uniprot ID:Q9Y6P5
Uniprot	Q9Y6P5
GeneID	27244;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SESN1

Note: This product is for in vitro research use only