

## Splicing factor, proline- and glutamine-rich

Catalog No: #AP87911



Package Size: #AP87911-1 50ug #AP87911-2 100ug #AP87911-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Splicing factor, proline- and glutamine-rich
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNA-binding p52/p100 complex, 100 kDa subunit, Polypyrimidine tract-binding protein-associated-splicing factor, PSF, Psf
Accession No.	Uniprot ID:Q8VIJ6
Uniprot	Q8VIJ6
GeneID	71514;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Sfpq

Note: This product is for in vitro research use only