

Myeloblastin

Catalog No: #AP87915



Package Size: #AP87915-1 50ug #AP87915-2 100ug #AP87915-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myeloblastin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proteinase 3,PR-3
Accession No.	Uniprot ID:Q61096
Uniprot	Q61096
GeneID	19152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prtn3

Note: This product is for in vitro research use only