

## RuvB-like 1

Catalog No: #AP87918



Package Size: #AP87918-1 50ug #AP87918-2 100ug #AP87918-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RuvB-like 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	49 kDa TATA box-binding protein-interacting protein, 49 kDa TBP-interacting protein, DNA helicase p50, Pontin 52, TIP49a, Tip49, Tip49a
Accession No.	Uniprot ID:P60123
Uniprot	P60123
GeneID	65137;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene\_full\_name: Ruvb1

Note: This product is for in vitro research use only