

NAD-dependent deacetylase sirtuin-1

Catalog No: #AP87921



Package Size: #AP87921-1 50ug #AP87921-2 100ug #AP87921-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NAD-dependent deacetylase sirtuin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hSIRT1,SIR2-like protein 1,hSIR2,SIR2L1,Regulatory protein SIR2 homolog 1
Accession No.	Uniprot ID:Q96EB6
Uniprot	Q96EB6
GeneID	23411;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SIRT1

Note: This product is for in vitro research use only