

Mothers against decapentaplegic homolog 2

Catalog No: #AP87925



Package Size: #AP87925-1 50ug #AP87925-2 100ug #AP87925-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mothers against decapentaplegic homolog 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MAD homolog 2, Mad-related protein 2, SMAD family member 2, SMAD 2, Madh2
Accession No.	Uniprot ID:O70436
Uniprot	O70436
GeneID	29357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Smad2

Note: This product is for in vitro research use only