

Signal transducer and activator of transcription 3

Catalog No: #AP87928



Package Size: #AP87928-1 50ug #AP87928-2 100ug #AP87928-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acute-phase response factor,APRF
Accession No.	Uniprot ID:P40763
Uniprot	P40763
GeneID	6774;
SDS-PAGE MW	22
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-200aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STAT3

Note: This product is for in vitro research use only