

Serpin B10

Catalog No: #AP87931



Package Size: #AP87931-1 50ug #AP87931-2 100ug #AP87931-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serpin B10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Bomapin,Peptidase inhibitor 10,PI-10,PI10
Accession No.	Uniprot ID:P48595
Uniprot	P48595
GeneID	5273;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPINB10

Note: This product is for in vitro research use only