

## Kit ligand

Catalog No: #AP87932



Package Size: #AP87932-1 50ug #AP87932-2 100ug #AP87932-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kit ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Mast cell growth factor, MGF, Stem cell factor, SCF, c-Kit ligand, SCF
Accession No.	Uniprot ID:Q28132
Uniprot	Q28132
GeneID	281885;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: KITLG

Note: This product is for in vitro research use only