

Stress-induced-phosphoprotein 1

Catalog No: #AP87933



Package Size: #AP87933-1 50ug #AP87933-2 100ug #AP87933-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stress-induced-phosphoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ST11,Hsc70/Hsp90-organizing protein,Hop
Accession No.	Uniprot ID:Q3ZBZ8
Uniprot	Q3ZBZ8
GeneID	617109;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STIP1

Note: This product is for in vitro research use only